Spectroscopic Determination of Excitation Temperature and Electron Density in Premixed Laminar Flame

S.S. Hamed

Physics Department Faculty of Girls, Ain Shams University, Cairo Egypt

Spectroscopic diagnostic technique are applied for the determination of excitation temperature T_{exc} and electron density in Air - Acetylene premixed laminar flame. A group of FeI spectral lines has been used to estimate T_{exc} using Boltzmann plot. Electron number density was determined from atom ion line pairs for Cd and Zn using Saha - Boltzmann equation after substitution the obtained value of excitation temperature. The obtained result of temperature show that the populations of energy levels follow Boltzmann's distribution law.

1. Introduction:

The atomization processes taking place in flames used as atom reservoir in analytical atomic absorption spectroscopy is considered as essential parameter that determines the absorption signal. The conditions with the flames depend mainly on the temperature and electron number density prevailing in the flame. These physical quantities describing the flame condition depend on some experimental parameters namely; fuel oxidant ratio, flow rates and design of the burner.

In flames burning at atmospheric pressure the primary reaction zone is extremely thin. Thus the residence time of the gases in this region is too short for the reaction between fuel and oxidant to reach completion or for thermodynamic equilibrium to be established. The region of interconal zone is much broader than the reaction zone; in such region the equilibrium is usually established via recombination of radicals and dissociation of the combustion products. This region of thermal equilibrium begins almost outside the reaction zone. The nature of the interconal zone for particular fuel / oxidant combinations varies considerably with the mixture strength. Consequently, this region is of interest for studies of the kinetics of flame processes, which depend mainly on excitation temperature in the flame. For the determination of the excitation temperature in flames sodium line reversal method is usually employed. [1,2,3] In this method sodium is added to the flame that is viewed through a spectroscope against a bright background continuum source. If the latter is operated to be hotter or colder than the flame, the sodium D lines will appear respectively either as dark lines in absorption or as bright light lines standing out against the the background. A reversal point is obtained when both source and flame are at the same temperature. This method requires special apparatus and are frequently time – consuming.

Another widely used technique is the intensity measurement of spectral lines emitted by the flame, which, in case of local thermal equilibrium, is given by :

$$I_{ij} = \frac{1}{4\pi} \frac{hc}{\lambda} \frac{N_o}{U(T)} g_i A_{ij} e^{-\frac{E_i}{KT}}$$
(1)

where N_{o} is the total number density of atoms in the ground state.

 g_i the statistical weight of the upper level.

 E_i is the energy of the upper level

U(T) is the partition function.

when oscillator strength f_{ij} is introduced instead of transition probabilities A_{ij} , the equation takes the form

$$I_{ij} = h \frac{e^2}{2\lambda^3} \frac{N_o}{U(T)} \frac{g_i f_{ij}}{\varepsilon_o m_e} e^{-\frac{E_i}{KT_{exc}}}$$
(2)

where ε_{a} is the permeability of free space.

 m_{e} the electron mass.

e the electron charge.

According to the forgoing equation T_{exc} can be determined using only one spectral line. However, this requires the determination of the absolute intensity, accurate values of transition probability, partition function and the knowledge of total number of atoms N_0 . To avoid these difficult T_{exc} can be determined from the intensity ratio of two or more spectral lines belonging to the same ionization stage. In this case, the partition function, and the number density N_o cancel, in addition there is no need to measure absolute intensity. In order to obtain accurate temperature values a number of spectral lines can be used. From the foregoing equation it can be seen that

$$log \frac{I\lambda^3}{g_f} = \text{Const.} - \frac{\text{E}_{i}}{\text{KT}_{\text{exc}}}$$

Thus the plot of $\log \frac{I\lambda^3}{gf}$ or $\log \frac{I\lambda}{gA}$ versus the energy of the upper

level E_i yield a straight line, if the population of the energy states obey Boltzmann distribution law. The straight lines is called Boltzmann plot. Its slope equals $-\frac{1}{KT_{exc}} = -\frac{0.625}{T_{exc}}$ if E_i is taken in cm⁻¹ or -5040/T if E_i is in eV.

The electron density, on the other hand can be determined by apply Saha equation from which the ratio of the intensity of an atom line I_a to that of an ion I_i belonging to the next ionization stage of the element is given by

$$\frac{I_{i}}{I_{a}} = \frac{1}{n_{e}} 2\left(\frac{2\pi nk}{h^{2}}\right)^{3/2} \frac{g_{i}A_{i}\lambda_{a}}{g_{a}A_{a}\lambda_{i}} T_{ion}^{3/2} \exp\frac{-E_{ion} + \Delta E_{ion}}{KT_{ion}} \exp\frac{E_{a} - E_{i}}{KT_{exc}}$$
(3)

From which the electron density (n_e) can be determined.

In the present work spectroscopic diagnostic technique are applied for the determination of excitation temperature and electron density in premixed laminar flame. The excitation temperature is measured using group of Fe I spectral lines [4] while the electron density is determined from the intensity ratio of atom / ion line pairs of Cd and Zn.[5]

2. Experimental:

2.1. Experimental Set up:

Perkin-Elmer atomic absorption spectrometer model (2280) equipped with burner nebulizer for an air - acetylene flame having single slot 100 mm was used. The dispersive element of the monochromator is plane grating having 1800 line/mm blazed at 255 nm. The monochromator has a focal length equal 267 mm equipped with a photomultipler and a digital read out system.

2.2. Spectral Lines for Excitation Temperature:

Table (1) gives the wavelengths and the atomic constants of 7 Fe I spectral lines. For These lines at least 12 sources for the oscillator strength values are available. For the determination of T_{exc} of the present flame cell these sets of oscillator strengths were applied to construct the Boltzmann plot. For this purpose a computer linear regression program (Segma plot) was used, that gives the slope of the straight line Boltzmann plot as well as the error in the slope ΔS that originates from the scatter of the experimentally determined line intensities. Accordingly both T_{exc} and the error ΔT_{exc} could be determined for each set of gf values, the most appropriat set of gf values was chosen that yields the smallest values of ΔT_{exc} .

Table (1): Fe I emission Line data.

		Log gf											
λ (nm) E (cm ⁻¹)		СТ	CB	MS	С	R _{etal}	K	AA	J	VS	BK	HP	BH
371.9	26875	-0.32	-0.29	-0.54	-0.54	-0.17	-0.54	-0.51	-0.27	-0.20	-0.43	-0.43	-0.45
372.3	27560	-1.02	-1.08	-1.28	-1.33	-1.02	-1.18	-1.22	-0.92	-1.03	-1.28	-1.22	-1.28
273.5	33695	0.57	0.62	0.35	0.40	0.57	0.41	0.38	0.67	0.52	0.31	0.30	0.23
273.7	27167	-0.50	-0.49	-0.66	-0.67	-0.31	-0.62	-0.64	-0.35	-0.35	-0.57	-0.57	-0.60
374.8	27560	-0.87	-0.83	-1.04	-1.08	-0.76	-1.00	-1.04	-0.72	-0.79	-1.01	-0.98	-1.06
374.9	34040	0.44	0.43	0.20	0.24	0.41	0.30	0.29	0.57	0.39	0.17	0.18	0.09
375.8	34329	0.31	0.34	0.03	0.07	0.22	0.18	0.31	0.45	0.26	0.00	0.00	-0.03

CT: Corliss & Tech [6]AA: Allen & Asaad [12]CB: Corliss & Bozman [7]J: Jurgens [13]MS: Margoshes & Scribner [8]VS: Valters & Startsev [14]C : Crosswhite [9]BK: Bridges & Kornblith [15]Retal: Reif & Fassel and Knisely[10]HP: Huber & Parkinson [16]K: King [11]BH: Banfield & Huber [17]

3. Spectral Lines for Determination of Electron Density:

For the determination of the electron number density in the flame under investigation, atom /ion line pairs of elements that have measured values of transition probabilities were selected. In the present work atom / ion line pairs of Cd and Zn will available. There E and gA values are given in Table (2). For these line pairs transition probabilities determined by two outhers are given in Table (2).

Species λ (nm)		E (cm ⁻¹)	$g_a A_a \lambda_i / g_i A_i \lambda_a$ Ratios			
			AS	BS		
Cd I	228.802	43692	2.6575	3.3949		
Cd II	226.502	44136	2.0375	3.3747		
Zn I	213.86	46745	3.0290	2.5203		
Zn II	206.19	48481				

 Table (2): Cd & Zn Data for electron density measurement.

AS: Andersen & Sorensen [18]

BS: Bauman & Smith [19]

4. Results and Discussion:

4.1. Optimal Condition:

The flame under investigation is used as an atom cell for atomic absorption spectroscopy measurements [20, 21]. For this reason it is worth to carry out that the spectroscopic diagnostic measurements under the optimal conditions of the flame as an atom reservoir. The parameters that control the flame conditions are the flow rates of fuel (acetylene), oxidant (air) and height of burner.

Optimal values of these parameters are considered those that yield the high value of absorption. It was observed that the optimal values for the flow rate of air and acetylene are 9.35L/min and 0.7L/min respectively; the values of the vertical and horizontal burner height are 17.5 mm and 20mm respectively.[22]

4.2. Atom Excitation Temperature:

Solution containing 200 µg/ml Fe was aspirated into the flame set at the previous obtained optimal conditions for the determination of T_{exc} . Table (3) gives the results obtained for the excitation temperature T_{exc} from the measuring intensity of the selected Fe spectral line after the substitution with all available sets of transition probabilities. It can be seen that two values of calculated T_{exc} may have nearly the same values, while the correction ΔT may differ significantly according to the set of transition probability. Figures (1&2) show an example for the Boltzmann plot obtained using the *gf* values according to Jurgens [13] and King [11]

Reference of <i>gf</i> values	Temperature (T)	ΔT_{exc}
СТ	2485.89	68.94
СВ	2607.54	98.03
MS	2555.20	66.07
С	2457.28	83.43
R _{etal}	2743.20	125.49
К	2719.52	473.46
AA	2392.54	105.17
J	2465.09	50.77
VS	2707.27	119.60
ВК	2710.14	122.48
HP	2745.45	97.70
BH	2776.09	128.99

Table (3): Results obtained for flame temperature measurements




Fig. (1): Boltzmann plot according to King values of gf.




Fig. (2): Boltzmann plot according to Jurgens values of gf.

5. Electron Density Measurement:

For the determination of electron density in the flame, standard solution containing 200µg/ml for Cd and Zn was nebulized at the previous optimal condition. The intensity of atom (I) and ion (II) lines for the two elements chosen was measured. The calculation of the electron density was carried out under the assumption that the flame is in local thermodynamic equilibrium (LTE). So the substitution for the ionization temperature (T_{ion}) by the same value obtained for the excitation temperature (T_{exc}), according to Jurgens g_f values, in equ.(3) the electron density for the two sets of the available transition probabilities. It can be seen from the table that the results are in good agreement with each other.

Table (5): Results obtained	for electron density	y measurement
-----------------------------	----------------------	---------------

Species	Reference of gA values	Electron density (cm ⁻³)
Cd	AS	4.770×10^{13}
Zn	AS	0.862×10^{13}
Cd	BS	3.750×10^{13}
Zn	CD S	1.039×10^{13}

Conclusions:

- From the results obtained in the present work it can be seen that straight lines obtained for the Boltzmann plots for *Fe* I lines show that the population distribution over the energy levels follow Boltzmann's law.
- In the Air- C_2H_2 flame the excitation is thermal with a temperature around 2400K⁰.
- The obtained excitation temperature was in good agreement with that obtained experimentally by continuum atomic absorption method. [23]
- The calculated electron densities for Cd from the two sets of gA used are in good agreement with each other as well as from Zn.

References:

- 1. A.G. Gaydon and H.G. Wolfhard, "*Flames, Their Structure, Radiation and Temperature*", Chapman and Hall, London, 3 rd ed., (1970).
- 2. R.H. Tourin, "Spectroscopic Gas Temperature Measurements", Elsevier, Amsterdam, (1966).
- **3.** G.F. Kirkbright, M. Sargent and S. Vetter, *Spectrochim. Acta*, **25B**, 465 (1970).
- 4. M.A.Eid, El Bialy, M.H. Abdallan, H.H.Higazy and S.S. Hamed, *Canadian J. Spectroscopy*, **35**, (1), (1990).
- 5. D.J. Kalnicky, V.A. Fassel and R.N. Kniselely, Applied Spec. 31 (1977).
- 6. C.H. Corliss and J.L. Tech "*Transition Probabilities for 3228 lines of Fe I*," Nat. Bur. Stand., Monograph 108 (1968).
- 7. C.H.Corliss and W.R. Bozman "Experimental Transition Probabilities for Spectral lines of Seventy Elements", Nat. Bur. Stand. Monograph 53 (1962).
- 8. M.Margoshes and B.F. Scribner, J.Res. Nat. Bur. Stand., 67A, 561 (1963).
- **9.** H.M. Crosswhite, "*The Spectrum of Fe I*", The Johns Hopkins University Press, Baltimore, (1958).
- **10.** I.G. Reif, V.A. Fassel and R.N. Kinseley, Reprints 16th Coll. Spectrosc. Int., Heidtberg, vol. 2, Hilger, London, (1970).
- 11. R.B. King and A.S. King, Astrophys. J. 87, 24 (1983).
- 12. C.W.Allen and A.S.Asaad, *Monthly Notices Roy. Aastron. Soc.* 117, 36 (1657).
- 13. G. Jurgens, Z. Physik, 138, 613 (1954).
- 14. A.K. Valters and G.P. Startsev, Opt. Spectrosc. (USSR), 17, 265 (1964).
- **15.** J.M. Bridges and R.L. Kornblith, *Astrophys. J.* **192**, 793 (1974).
- 16. M.C.E. Huber and W.H. Parkinson, Astrophys. J. 192, 793 (1974).
- 17. F.B. Banfield and M.C.E. Huber, Astrophys. J. 186, 335, (1973).

- **18.** T. Andersen, J. Desesquelles, K.A. Jensen, and G. Sorensen, J. Opt. Soc. Am. **60**, 1199 (1970).
- 19. S. R. Bauman and W.H. Smith, J. Opt. Soc. Am. 60, 345 (1970).
- **20.** R.Sinem Atgin, O.El-Agha, A. Zararsiz, A.Kocatas, H. Parlak, and G.Tuncel, *Spectrochimica Acta Part* **B55**, 1151 (2000).
- 21. U. Drivikli and L. Elci, Analyti. Chim. Acta, 452 (2), 231 (2002).
- **22.** A.B. El Bialy, S.S. Hamed, W.M. Moussa and R.K. Abd El-Hameed, *Egypt. J. Solids* (**28**), 151 (2005).
- 23. R. F. Browner and J. D. Winefordner, *Analytical Chemistry* 44 (2). 247 (1972).